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DE **DENAS®·Vertebra**
(s.67–100)

TABLE OF CONTENTS

1. Safety measures.....	38
2. Purpose of the device.....	41
3. General information on method, used in the device.....	42
4. Device description.....	43
5. Indications and contraindications for use.....	45
6. Conditions and order of treatment.....	50
7. Technical maintenance and order of replacement of batteries.....	56
8. Troubleshooting list.....	58
9. Specifications.....	60
10. Complete set.....	63
11. Manufacturer's warranty.....	64
12. Address of the manufacturer.....	65
Warranty maintenance form.....	101
Acceptance certificate.....	104

Thank you for purchasing device DENAS®-Vertebra.

We hope that this device will help you to maintain your GOOD health in many years to come.

In order to make use of the unit effective and safe, please, carefully read all sections of this manual.

1. SAFETY MEASURES

Pay attention to all information, marked with this sign. It is important for your safety and proper use of the device.

Information, contained in present operations manual, is important for your safety and proper use and maintenance of the device.

The device is safe for use, because it utilizes internal power source of low voltage, which is isolated from work part of the device (article of type B with body of type F).

The device must not be used for treatment of patients with implanted electronic devices (for example, pacemakers) and for treatment of patients who have individual electric current intolerance.

Use of the device in direct front projection of heart is prohibited.

Don't treat patient with any high-frequency electric device during stimulation; simultaneous use of the device and other electric equipment can cause burns and lead to possible damage of the device.

Work near short-wave and microwave equipment can bring to instability of output parameters of the device.

It is prohibited to use batteries, which are different from specified by the manufacturer. Use of inappropriate batteries can lead to device malfunction, and also can result in electric shock.

It is prohibited to deform module of electrostimulation.

The device contains fragile components. Protect it from shocks.

Protect the device from influence of direct sunbeams under high ($>25^{\circ}\text{C}$) air temperatures. It can lead to malfunction of the device indicator or destroy batteries.

The device is not waterproof. Protect it from ingress of moisture.


All works on maintenance and repair of the device must be executed by qualified specialists of the manufacturer.


Attention! Prevent the device or its parts (control panel, utility power source) from falling on the ground.


Attention! If the patient is running through a course of other physical therapy, use of the device is possible only under approval of attending physician.

Transportation conditions: temperature from -50°C to $+50^{\circ}\text{C}$, relative air humidity from 30 to 93%, atmospheric pressure from 70 to 106 kPa (525 to 795 mm Hg)

Storage conditions: temperature from -50°C to $+40^{\circ}\text{C}$, relative air humidity from 30 to 93%, atmospheric pressure from 70 to 106 kPa (525 to 795 mm Hg)

Operational conditions: temperature from +10 °C to +35 °C, relative air humidity from 30 to 93%, atmospheric pressure from 70 to 106 kPa (525 to 795 mm Hg).

Attention! If the device has been stored at the temperature below 10 °C, keep it in normal climate conditions for no less than 2 hours - before use.

Utilization: All packaging materials are not environmentally harmful, they may be used repeatedly.

Separate collection of electrical and electronic equipment.

Attention! The device contains valuable materials, which can be used repeatedly after utilization with consideration of requirements of environmental protection. They shall be delivered to specially intended for this purpose places (consult with corresponding services in your district) for collection and processing.

2. PURPOSE OF THE DEVICE

The device DENAS-Vertebra is intended for conduction of dynamic electrostimulation and electromassage by way of impact on different zones and areas of the back, and also for pain relief and correction of functional disturbances under diseases of different body organs.

The device is intended for use in hospitals and in home conditions in accordance with the instructions of a physician.

Compliance with standards: This medical device is CE marked according to the Directive 93/42 /EEC on medical equipment.

3. GENERAL INFORMATION ON METHOD, USED IN THE DEVICE

One of the most promising areas of modern medicine is use of physiotherapy for treatment and rehabilitation. Currently, there goes major distribution of complex technologies, such as neurostimulation (dynamic electroneurostimulation). The method consists in the therapeutic effects of short asymmetrical bipolar current pulses, shape of which varies depending on the electrical resistance of tissues in zone under electrodes. Safe, simple and with minimal contraindications, method DENS is used in the treatment of a wide range of diseases. At present time DENS has been successfully applied in many fields of medicine; extensive data on the benefits of neurostimulation in treatment of various diseases of adults and children has been received. In addition, specially selected modes of stimulation of the device allow simulating massage effect, due to muscle contraction. The effects and feelings during electromassage are comparable with manual massage, and often surpass it due to deeper treatment of muscles.

4. DEVICE DESCRIPTION


Fig. 1. Assembled device. Top view


Fig. 2. Module of electrostimulation (bottom view of the part, which contains electronic unit)


Fig. 3. Keyboard of control panel


Fig. 4. Designation of symbols on the device indicator

5. INDICATIONS AND CONTRAINDICATIONS FOR USE

EN

General purpose of the device

The device can be used in following cases:

- as the primary method of treatment under drug intolerance and contraindications for application of other methods;
- as a component of complex treatment in order to enhance the effect of primary drug, homeopathic, manual therapy as well as mental therapy and other methods of treatment;
- as a symptomatic treatment under different diseases and syndromes.

Attention! First of all, the only sign of a serious illness often is the appearance of sudden pain in any part of your body. Therefore, if episodes of pain appear for the first time and repeat, and intensity of the pain increases, urgent consultation with a doctor is required.

Attention! The use of the device does not exclude need for other methods of treatment, prescribed by your doctor.

Attention! Impact, which is formed by the device, is not a therapeutic massage.

Indications for use:

- acute and chronic pain syndromes under diseases of spine, joints and muscles;
- pain relief and correction of functional disorders under diseases of internal organs;
- rehabilitation after previous diseases, surgical interventions, traumas;
- enhancement of adaptive capacity of organism under physical and mental strain.

Contraindications:

- idiosyncrasy;
- presence of an implanted pacemaker;
- status epilepticus;
- neoplasms of any etiology and distribution;
- vein thrombosis;
- acute agnogenic febrility;
- state of acute mental excitation of alcohol or drug origin.

Attention! There have been no controlled researches on impact of the device on pregnant women, so the device should be used only if it is absolutely necessary, after consultation with the attending physician.

Attention! In the case of presence of above contraindications, possibility of use of the device should be discussed with a doctor.

Attention! Use of the device in direct front projection of heart is prohibited.

Attention! It is prohibited to use the device in presence of solution of continuity of skin (wounds, skin diseases) in zone of impact.

Programs of the device

Pro- gram	Main zone of impact	Area of application	Period of program operation (min.)
1	2	3	4
A	Zone of the back at full length of the spine	<p>1. Treatment, rehabilitation, complementary prevention of chronic diseases.</p> <p>2. Prevention and treatment of stress.</p> <p>3. Enhancement of adaptive capacity of organism under intense physical and mental work, physical and mental fatigue, chronic fatigue syndrome, difficulties under waking in the morning and afternoon sleepiness, sleep disorder in the evening and insomnia, hyperirritability, sexual function disorder.</p>	23 ± 1
B	Zone of the back at the level of cervical spine, thoracic spine or lumbar spine - optional	<p>1. Severe topalgia in the back, associated mainly with rachiopathy and damage of back muscles.</p> <p>2. Spastic pain under diseases of internal organs.</p> <p>3. Rehabilitation in the postoperative period</p>	23 ± 1

1	2	3	4
C	Zone of the back at the level of cervical spine, thoracic spine or lumbar spine - optional	Moderate and mild pain in the back, associated with spine injuries. Functional disorders, associated with overexertion of back muscles: fatigue, "choking" of muscles, etc. Chronic pain syndromes under diseases of spine and internal organs. Correction of functional disorders under internal diseases. Treatment, rehabilitation at all stages, prevention of long-term illness of the spine.	26 ± 1
D	Impact on the selected zone of the back with a small area (applicator mode)	Emergency care under acute back pain, associated with rachopathy and / or diseases of spinal nerves. Improvement of effectiveness of other programs	10 ± 1

Attention! If a patient suffers from urolithiasis, before use of the device it is advised to consult with a specialist (urologist).

Attention! Insofar if in the list of indications there is no data on some disease, then before use of the device, it is recommended to consult with a specialist.

Attention! It is not recommended to use the device for more than 40-50 minutes in one session. Exceeding time-limit is allowed only under presence of doctor's pre-trial detention.

Effectiveness of treatment

Under use of the device, it is necessary to consider that the most intense therapeutic effect is a result of a course of treatment. Its duration for some diseases should be 6-8 sessions, for other diseases - 8-12 sessions, rarely 14-20 sessions. Total therapeutic effect of separate procedures also ensures a long after-effect of the treatment, which continues on after its completion, and in some cases, long-term results are even more significant than direct results. However, keep in mind that prolonged continuous use of a physical factor increases adaptation to it and may reduce the effectiveness of a treatment. In addition, the probabilistic nature of physiological processes in the body leads to lack of effectiveness for a small number (5%) of people or lack of favorable reactions from use of therapeutic factors.

Attention! *In the case of breakdown of health during treatment, you should consult with a doctor.*

6. CONDITIONS AND ORDER OF TREATMENT

6.1. General Provisions

For applying of treatment with the help of the device, there is no need to observe any special conditions. The device should be placed in a warm room on a hard flat surface (medical couch, floor), length and width of which should correspond to that of the patient. After the procedure, it is recommended to relax for 10-15 minutes.

Patient should be located on the module of electrical stimulation in such a way that majority of the electrodes along the length of patient contact with the skin surface. Otherwise (under absence of skin contact or under contact with less than 60% of electrodes) program may not start or stop at one of its stages. Most of programs of the device, regardless of selected zone, have one or more phases of stimulation along whole surface of the back. In case when during this phase part of electrode groups will not contact with the skin, the device will abort program due to impossibility of its further implementation.

When properly positioned on the module of electrical stimulation, upper paired electrodes are in the zone of neck just below the hairline, and central electrodes are in contact with your skin in projection of the spine. In case of inconvenience, associated with hyperextension of cervical spine, small pillow, folded towel or similar object can be placed under the head.

Automation of the process of treatment provides ease and convenience of procedures conduction. During the session, the device switches electrodes and stimulation frequency in order, specified by the program. Combination of certain switching speeds and frequency of impact secure simulation of massage effect, which is most fully manifested under use of comfortable and maximum power levels.

Attention! Programs of the device consist of several phases, which subsequently change each other. Each phase has different parameters of stimulation, including repetition rate of impulses. Under changing of phases, subjective feeling of change of impact power level is possible. In case of appearance of sense of discomfort, adjust the power, using the appropriate buttons on control panel.

6.2. Determination of intensity (power) of impact

The level of energy treatment with impulse current is determined individually based on subjective sensations of the patient. Intensity of electrostimulation is conditionally divided into three levels.

First, *minimal level* — the patient does not have any subjective sensations or has subtle vibrations in area of impact under electrodes. This level is used when treatment shall not be intensive –for patients with diseases of cardiovascular system (hypertension, hypotension, a syndrome of vegetative-vascular dystonia). For patients who are experiencing frequent headaches, dizziness, treatment should also be carried out with minimal power.

Second, *comfortable level* — patient feels light pricking, vibration or light burning without pain. It is used under treatment of weak pain and aches of medium intensity. This is the most frequently used power level.

The third, *maximal level* — the patient feels expressed painful pricking or burning. Such intensity can be accompanied by involuntary contraction of muscles close to electrodes. It is used under expressed pain syndrome.

Attention! Control of power level is executed subjectively, considering sensations of patient. You shall not exceed pain threshold.

Attention! For patients with hypertension, hypotension, a syndrome of vegetative-vascular dystonia, as well as patients who have frequent headaches, dizziness, treatment should be carried out with minimal power level.

Attention! Power of electrostimulation can be increased or reduced during the treatment course depending on the extent of the patient's sensitivity

6.3. Preparation for treatment

Open consumer's container and get the device. If there is no battery in the battery compartment, set batteries. Under setting of batteries, observe polarity. Before use of the device, make sure that on the body of control panel, as well as on module of electrostimulation, there are no traces of falls, shocks, LCD display is not broken, connecting cable is not damaged. Make sure that there are no traces of foreign substances and traces of damage on electrodes.

Attention! It is prohibited to use the machine if it is damaged.


Attention! If the device has been stored in cold environment or in high humidity conditions, keep it under normal room temperature for no less than 2 hours - before use

Attention! Treat the device electrodes with a standard disinfection means (such as 3% solution of hydrogen peroxide) after each procedure).

6.4. Order of conduction of treatment session

6.4.1. Lay-out of sheet

Get from the packaging one sheet, specifically designed for the device. Place it on electrode surface of the module of electrostimulation in such a way that holes of the sheet fully match location of electrode groups. Place the sheet below level of electrode groups, so that it will be drawn tight and securely fastened.

6.4.2. Switching on the device

Prepare the device for conduction of procedure. Press button on control panel. The device will switch on and switch to standby mode.

6.4.3. Selection of impact program

Upon switching on the device, program A is used by default. To select required impact program press button

④. Program are located in the menu in following order:
program A → program B → program C → program D →
program A...

6.4.4. Selection of impact zone

For programs, which are associated with necessity of selection of impact zone (programs B, C and D), use button ⑤. If indicated on the display zone does not correspond to area of wanted stimulation, shift stimulation zone up and down, depending on your choice.

6.4.5. Start of program and setting of power level

Comfortably position yourself on the module. Under contact with skin of more than 60% of electrodes in selected zone, there will be shown a symbol ⑥ – on presence of contact of the device with skin

Start of the program is executed by button ⑦ under presence of contact of electrodes with skin. To select power level, use buttons ⑧ and ⑨:

- single depression changes power level for one point;
- holding of a button for more than 1 second will lead to consistent change of impact level. After selection of required level, release the button. Countdown of impact time will begin.

After the end of the program, device will emit sound signal and stimulation will end.

Also the device emits signal, when there is no contact of electrodes with skin for more than 5 seconds (less than 60% of contact group of selected zone).

6.4.6. Volume setting

For convenience, you can set volume of signaling. The device can work without sound or under two levels of volume of sound signals. Volume can be adjusted by pressing two buttons ⑩ and ⑪ or ⑫ or ⑬.

on and increase of volume level simultaneous pressing of combination of buttons  and , and for decreasing of volume or switching the sound off — simultaneous pressing of combination of buttons  and . Volume is displayed on the device indicator, image is absent; ;).

6.4.7. Switching off the device

Press and hold button  until the device will be switched off.

6.4.8. After the end of session, it is recommended to relax for 10-15 minutes.

6.4.9. Completion of procedure

For switched off device, remove napkin and treat electrode surface of unit of electrostimulation with standard disinfectant (e.g. 3% solution of hydrogen peroxide). In case of serious dirt accumulation on the module, it is allowed to treat the device with lint-free cloth, moistened with neutral and non-abrasive detergent.

Attention! Lack of adequate cleaning of elements can lead to risk of allergic reactions appearance or cross colonization - if the device is used by several people.

Attention! In case of idiosyncrasy, if during a session, there is deterioration of health or after a session you feel skin irritation, it is necessary to stop treatment and consult a doctor immediately.

Attention! Under appearance of any allergic reactions in zone of contact of the device with skin, you should stop using the device and consult a doctor immediately.

7. TECHNICAL MAINTENANCE AND ORDER OF REPLACEMENT OF BATTERIES

Technical maintenance of the device should include following procedures:

7.1. External examination of the product. You must make sure that there are no traces of blows and falls, because it can result in incorrect operation of the device.

7.2. Before using the device, and also after its use, it is necessary to clean module of electrostimulation. Lack of adequate treatment may cause allergic reactions or infection of the skin if the device is used by several people. To clean it, use standard disinfectants (e.g. 3% solution of hydrogen peroxide) and soft, lint-free cloth.


7.3. If you intend not to use the device for a long period of time (more than 3 days), remove the batteries from the battery compartment. This will prevent their early discharge.

7.4. Under flashing of a symbol of battery status on LCD display, you should replace batteries. To do this:


- take off the cover of battery compartment, as shown on fig. 5;
- take off used batteries;
- put new batteries, observing polarity;
- fix cover of battery compartment back into its place.


Attention! Use only high-quality batteries! If you use defective batteries, there is a danger of spontaneous discharge of batteries contents, as a result, the device may become unserviceable (not covered by manufacturer's warranty), and there is a risk of chemical burn.

1. Take-off the cover of battery compartment and get used batteries.

2. Put new batteries, observing polarity.

3. Fix cover of battery compartment back into its place.

Fig. 5. Order of replacement of batteries

8. TROUBLESHOOTING LIST

In this table conditions of the device are listed – which can be treated as malfunction and corrected independently. In case of other malfunctions, contact with representative of the manufacturer (see section 11 «Manufacturer's warranty»), don't try to correct them by yourself.

Trouble	Possible reason	Method of elimination
1	2	3
The device does not switch on when you press the button	There are no batteries	Set new batteries (see section «Technical maintenance and order of replacement of batteries»)
	Voltage of batteries is less than 2,1 V	
Under switching on, the device emits sound signals and switches off automatically	Voltage of batteries is less than 2,1 V Period of time between switching on of the device and its subsequent switching on is less than 3 seconds	Change batteries (see section «Technical maintenance and order of replacement of batteries») Wait for 3 seconds after switching off the device and switch it on again
The device switches off or symbol of battery flashes and the device emits discontinuous signal	Voltage of batteries is less than 2,1 V	Set new batteries (see section «Technical maintenance and order of replacement of batteries»)

1	2	3
The device does not pass to mode of electrostimulation under contact of built-in electrodes with skin. There is no symbol of contact of the deivce with skin on LCD-indicator	Insufficient area of contact of electrode with skin	Make sure that no less than 60% of total number of electrodes in each zone of the device are tightly applied to the skin. Take off your underwear
During work, the device emits sound signal, shuts down the program and exits to menu		
The device does not emit sound signals	Sound is switched off	Switch the sound on according to p.6.4.6
«Quick consumption» of batteries	Off-grade batteries	Use high-quality batteries (alkaline batteries are recommended) or accumulators of corresponding unitsize with voltage of no more than 1,5 V

9. SPECIFICATIONS

9.1. Maximum current consumption of the device under supply voltage of $3,0 \pm 0,15$ V:

- in mode of electrostimulation – no more than 300 mA;
- in switched off mode – no more than 100 uA.

9.2. Weight of the device without batteries:

- electrostimulation module – no more than 5,000 g;
- control panel – no more than 180 g.

9.3. External dimensions of the device are no more than:

- a) control panel - 140x55x28 mm;
- b) electrostimulation module- 900x365x70 mm.

9.4. The unit automatically shuts off after three minutes since last pressing of any button except button , when programs of the device are not used.

9.5. Characteristics of the device:

9.5.1. Program A

Phase	1	2	3	4	5	6	7	8
Frequency (Hz)	60	77	140	140	200	20	77	140; 77; 60; 20
Number of impulses per electrode	30	40	70	30	3	30	60	5
Phase length	3 min	Total time - 2 minutes. Every frequency works 30 seconds						

Program B

Phase	1	2	3	4	5	6
Frequency (Hz)	200	77	60	77	20	140; 77; 60; 20
Number of impulses per electrode	1	1	30	40	30	5
Phase length	10 min	5 min	2 min	2 min	2 min	Total time - 2 minutes. Every frequency works 30 seconds


Program C

Phase	1	2	3	4	5	6	7
Frequency (Hz)	60	77	140	140	20	77	140; 77; 60; 20
Number of impulses per electrode	30	40	70	30	30	60	5
Phase length	5 min	5 min	5 min	3 min	3 min	3 min	Total time - 2 minutes. Every frequency works 30 seconds

Program D

Phase	1	2	3
Frequency (Hz)	140	77	20
Number of impulses per electrode	1	1	1
Phase length	5 min	3 min	2

9.5.2. Parameters of separate impulses of impact:

Pulse form under connected load Z


Under minimum power level:

- duration of positive part of impulse (T1f) - 25 ± 10 us;
- amplitude of positive part of impulse (A1f) no more than 30 V;
- impulse swing (Ap-p) no more than 70 V.

Under maximum power level:

- duration of positive part of impulse (T1f) - 250 ± 70 us;
- amplitude of positive part of impulse (A1f) - 30 ± 10 V;
- impulse swing (Ap-p) - 190 ± 70 V.

Parameters of load Z.

R1 — $56 \text{ kOhm} \pm 10\% - 0,25 \text{ W}$

R2 — $470 \text{ Ohm} \pm 10\% - 0,25 \text{ W}$

C1 — K10 - 47a - H30 - 6200

pF $\pm 20\%$ - 250 V

10. COMPLET SET

Below you can see list of delivery for the device:

Item	Quantity, pcs.
Device DENAS-Vertebra	1
Operations manual	1
Safety sheet	10^2
Battery of type LR6/AA ¹	2^2
Consumer's container	1

¹ Manufacturer's warranty is not applied for these accessories.

² Quantity can be increased.

11. MANUFACTURER'S WARRANTY

11.1. Operating time of the device — 5 years.
Under observance of operation regulations, actual operating time can significantly exceed official.

11.2. Warranty period for the device — 24 months from the date of sale.

11.3. The seller (manufacturer) or organization, which acts as a seller (manufacturer) on the basis of contract with it, is not responsible for the defaults should they occur after sale of the device as a result of:

- 1) failure on the part of the consumer to comply with the rules of transportation, storage, care and operation, provided by the present manual;
- 2) mechanical damages;
- 3) actions of the third party;
- 4) force-majeure.

11.4. Guarantee obligations do not apply to products with broken manufacturer's seal - which is located in batteries compartment.

11.5. In case of unit breakdown or malfunction within the warranty period, as well as in case of incomplete shipping, the owner must send the following documents to the manufacturer's address or manufacturers' representative: claim for repair (exchange) with name, address, telephone number; defects list with brief description of the malfunction, date and conditions of its appearance.

12. ADDRESS OF THE MANUFACTURER

Address of company-manufacturer:

LLC «RC ART»
Academika Postovskogo st., 15
Ekaterinburg, 620146, Russia
Telephone: +7 (343) 267-23-30
<http://www.denascorp.ru>
e-mail: corp@denascorp.ru

Official Representative in the European Union:

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